

Montage COLLECTION



shown in Oatmeal

Oatmeal



Brown Sugar



Ruby



Royal



Silt



Mahogany



Slip into the casual comfort of the all-new Montage Collection! Our new design features a plush, generously stuffed seam back, better seat support and a more casual armrest design for hours of relaxation. In four different sizes and six standard colors, the Montage collection is sure to become a household favorite!

Montage
542-JPT
For people 5'2"
and under

UltraComfort[™]
AMERICA
Positioning you for life.

PROUDLY

MADE IN AMERICA
with Domestic & Imported Materials

Montage COLLECTION

Montage UC 542-JPT

Model Number	UC542-JPT
Size	Petite
Weight Capacity	300 LB
Recline Positions	3
Chaise Seat	Yes
Pocketed Coil Seating System	No
Back Type	Seam
Head Rest Cover	Yes
Arm Covers	Yes
Battery Back Up	Yes
Overall Width	31"
Overall Height	39.5"
Width Between Arms	19"
Floor to Top of Seat	19"
Seat Depth	18.5"
Seat to Top of Back	26"
Seat to Top of Arm	6"
Seat Angle Raised	15°
Distance Required from Wall Reclined	17"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Frame, Recline Mechanism	Yes
Warranty on Electrical parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Warranty on Fabric and Foam, and Pocketed Coils	1 YR
Removable Back Option	No
Customizable *	Yes
Left Hand Control Option Available	Yes
Extra Pocket Option Available	Yes
4-Motor Heat & Massage Available	Yes
2-Motor Heat & Massage Available	No
Foot Rest Extension Option Available	No
Smart Tek System	Yes

* See Customize Chart for dimension limits
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness.
 All measurements taken on concrete floor using levels and metal rulers.



User-friendly feature of every UltraComfort lift and recline chair:

Each UltraComfort lift and recliner has a battery backup system in the electrical transformer which operates on two 9 volt batteries. In the event of a power failure, the battery system will raise the chair from the reclined to the seated position. The batteries should be replaced after one use.



UltraComfort
 AMERICA™

Positioning you for life.

UltraComfort
 Phone: 866-793-5100
 Fax: 866-255-9962
 web: www.ultracomfort.net



Literature is current at time of printing. UltraComfort reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68

Listing # E112417
 FNAE 80-214A

